

## DT5521HE

### 4 Channel 6 kV/20 $\mu$ A Desktop HV Power Supply Module (USB/Ethernet)



## Features



- 4 independent channels in a compact desktop case
- 6 kV / 20  $\mu$ A output range
- Available with positive, negative or mixed polarity
- SHV coaxial output connectors
- Low Ripple
- 100 pA Current monitor resolution (with x10 Imon-Zoom: 10 pA)
- Remote control (USB 2.0 and Ethernet)
- Software Tools for easy channel management

## Description

The CAEN **Mod.DT5521HE** is a desktop module housing **4** independent High Voltage Power Supply Channels **6 kV, 20  $\mu$ A**. The unit is available with either positive or negative output polarity; mixed version with 2 positive and 2 negative channels is also available. HV outputs are delivered through **SHV** connectors.

### SHV connector

Radial R317580 HV coaxial connector for Mod. DT5519HE

Consult our **connectors reference page** for technical information.

The HV output RAMP-UP and RAMP-DOWN rates may be selected independently for each channel in the **1 ÷ 500 V/s** range with **1 V/s** steps. The module features 100 pA I<sub>mon</sub> resolution. Features include current monitor resolution extended by 10x in a Lower range (selectable via software). Functional parameters can be programmed and monitored via **USB 2.0** or **Ethernet**. A complete set of software tools is available to control these units, from low level libraries to graphical application software.

Safety features allows the module to perform as a current generator and includes:

<b>Channels</b>	can be enabled or disabled through the Global Interlock logic.
<b>Overvoltage and Undervoltage warning</b>	when the output voltage differs from the programmed value ( $\pm 0.5\%$ of the set value, $\pm 3$ V).
<b>Overcurrent detection</b>	when a channel attempts to exceed the programmed current limit, it signaled to be in "overcurrent" and enter in a TRIP status. The output voltage is varied to keep the current below the programmed limit for a programmable TRIP time, then the channel is switched of. If TRIP is set to "constant current mode", the channel behaves like a current generator.
<b>Hardware VMAX, IMAX</b>	Programmable VMAX, IMAX protection limit.
<b>Safety Board Interlock</b>	common Interlock logic for channels enable/disable and individual inputs signal for channel Kill function.

## Technical Specifications

### Package

Desktop module housed in a 154x50x164 mm<sup>3</sup> (WxHxD) alloy box; weight: ~800g

### No. of Channels

4

### Output Voltage

0 ÷ 6 kV

### Polarity

Positive / Negative / Mixed depending on purchased version; common ground

### Max. Output Current

20 µA

### Voltage Set Resolution

100 mV

### Voltage Monitor Resolution

50mV

### Current Set Resolution

0.5 nA

### Current Monitor Resolution

100 pA (high range) / 10 pA (low range)

### Current Set Maximum Value

21 µA

### VMAX hardware

0 ÷ 6100 V selectable for each channel

### VMAX hardware resolution

25 V

### VMAX hardware accuracy

2% of FSR

### Ramp Up/Down

1 ÷ 500 Volt/sec, 1 Volt/sec step

## Trip

Max. time an "overcurrent" is allowed to last (seconds). A channel in "overcurrent" works as a current generator; output voltage varies to keep the output current lower than the programmed value. "Overcurrent" lasting more than set value (1 to 9999) causes the channel to "trip". Output voltage will drop to zero either at the Ramp-down rate or at the fastest available rate, depending on Power Down setting; in both cases the channel is put in the off state. If trip= INFINITE, "overcurrent" lasts indefinitely. TRIP range: 0 ÷ 999.9 s; 1000 s = Infinite. Step = 0.1 s

## Voltage Ripple

20 ÷ 1000 Hz: <7 mVpp typical; < 10 mVpp maximum  
1 ÷ 20000 kHz: <2 mVpp typical; < 5 mVpp maximum

## Vmon vs. Vout accuracy

±0.05% of read ±1 V

## Vset Vs. Vout accuracy

±0.05% of read ±1 V

## Imon vs. Iout accuracy

high range: ±2% of read ±5nA; low range: ±2% of read ±500pA

## Iset vs. Imon Accuracy

high range: ±2% of read ±5nA; low range: ±2% of read ±500pA

## Imon / Temperature coefficient

max 100 ppm/C°

## Maximum output power

180 mW per channel

## Ordering Options

Code	Description	
WDT5521HEMAA	DT5521HEM - 4 Channel 6 kV/20 $\mu$ A Desktop HV Power Supply (USB/Ethernet) - Mixed	RoHS
WDT5521HENAA	DT5521HEN - 4 Channel 6 kV/20 $\mu$ A Desktop HV Power Supply (USB/Ethernet) - Negative	RoHS
WDT5521HEPAA	DT5521HEP - 4 Channel 6 kV/20 $\mu$ A Desktop HV Power Supply (USB/Ethernet) - Positive	RoHS

## Accessories

### A1484

Inhibit - Kill Signal BNC Adapter for HV Power Supply Modules

---

### A148x



Inhibit - Kill Signal BNC Adapter for HV Power Supply Modules

---

### HV CABLES



High Voltage Cable Assemblies

---


### A1483

Inhibit - Kill Signal BNC Adapter for HV Power Supply Modules

---

## Related Software

### GECO2020



General Control Software for CAEN HV Power Supplies

### CAEN Toolbox



Multi-Functional Software Suite for the Upgrade of Front-end Boards, Bridges and Power Supplies

## Related Software Libraries

### CAEN HV Wrapper Library



Library for CAEN Power Supply Control

## Related Products

### LabVIEW Driver (PSM - Power Supply Modules)



LabVIEW Instrument Driver for Power Supply Modules

# Gallery



**This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit [www.caen.it](http://www.caen.it) © Caen S.p.A - 2024**

**CAEN S.p.A.**

Via Vetraia 11  
55049 - Viareggio  
Italy

**Phone +39.0584.388.398**

**Fax +39.0584.388.959**

**info@caen.it**

**www.caen.it**

**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt, Germany

**Phone +49 (0)212.254.4077**

**Mobile +49 (0)151.16.548.484**

**info@caen-de.com**

**www.caen-de.com**

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA

**Phone +1.718.981.0401**

**Fax +1.718.556.9185**

**info@caentechnologies.com**

**www.caentechnologies.com**

**CAENspa India Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri West  
Mumbai, Maharashtra 400053, India

**info@caen-india.in**

**www.caen-india.in**

