

N957

8k Multi-Channel Analyzer



Features



- 1 input Multichannel Analyzer
- Fully computer controlled MCA
- 8K ADC, fast conversion time with linearization enhancement circuit
- PHA acquisition mode
- Suitable for HPGe, Nai(Tl), CdTe and other detector technologies
- USB 2.0 communication interface
- Libraries, Demos (C and LabView) and Software tools for Windows and Linux

Description

The **Mod. N957** is a 8k **Multi-Channel** (MCA) with USB port, housed in a 1-unit wide standard NIM module. The multichannel analyzer performs the essential function of collecting the data and producing output, in the form of converted value of input peaks. The managed input pulses are those produced by a standard spectroscopy amplifier: Gaussian, semi-Gaussian or square waves, unipolar (positive) or bipolar, in a range from 0 to 10 V, with a rise time greater than 0.1 μ s. The trigger can be made "on signal" (Auto Gate mode) or "external" (External Gate mode). In the first case a discriminator, with a settable threshold, enables the conversion. In the second case, an external gate is fed to the module, via front panel GATE In connector. The input channel has one peak amplitude stretcher, the output of which is digitised by a 13 bit fast (0.8 μ s) ADC featuring a sliding scale technique, to improve the differential non-linearity. Converted waveforms are stored into a 64 KSamples buffer memory. The unit hosts a USB2.0 port (also compatible with USB 1.1), which permits a simple control and data-acquisition via PC.

Libraries for Windows and Linux will be provided as well. Future firmware upgrade is possible via USB.

Technical Specifications

Packaging

One unit NIM mechanics

No. of ADC channels

1

Resolution

13 bit (8192 channels - 8064 valid if sliding scale enabled)

ADC Conversion time

0.8 μ s

Dead Time

4.8 μ s

LSB

1.22 mV

Gate

Signal must occur prior to and must extend for at least 0.2- μ s after the peak (in External Gate mode)

Differential Non-Linearity

< 1%
from 5% to 95% of input FSR (500 mV \div 9.5 V)

Integral Non-Linearity

< 0.065%
from 5% to 95% of input FSR (500 mV \div 9.5 V)

Gain Instability

< +150 ppm/ $^{\circ}$ C

USB port

Compatible with USB 1.1 and USB 2.0;
3m maximum cable length (longer distance can be achieved with commercial off-the-shelf products)

Maximum transfer rate

30 Mbyte/s (USB 2.0 Bulk Transaction Protocol);
75 Kbytes/s (USB1.1)

I/O signals

NIM/TTL; selected via internal switch

Discriminator Threshold

Software programmable, 0 mV \div 500 mV range, 100 steps

Ordering Options

Code	Description
WN957XAAAAAA	N957 - 8k Multi-Channel Analyzer RoHS


Related Products

NIM8306



2 Slot Switching 750 W Mini Crate

NV8020A




7U CRATE VME/NIM 8 slot VME64 365W, 5 slot NIM 150W

NIM8303



5U 12 slot 300/600 W Crate

NIM8304



7U 12 slot smart fan unit Switching 2000 W Crate

NIM8305



2 Slot Switching 450 W Mini Crate

NIM8302



5U 10 slot 150 W Compact Crate

Gallery



This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it © Caen S.p.A - 2024

CAEN S.p.A.

Via Vetraia 11
55049 - Viareggio
Italy

Phone +39.0584.388.398

Fax +39.0584.388.959

info@caen.it

www.caen.it

CAEN GmbH

Brunnenweg 9
64331 Weiterstadt, Germany

Phone +49 (0)212.254.4077

Mobile +49 (0)151.16.548.484

info@caen-de.com

www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA

Phone +1.718.981.0401

Fax +1.718.556.9185

info@caentechnologies.com

www.caentechnologies.com

CAENspa India Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri West
Mumbai, Maharashtra 400053, India

info@caen-india.in

www.caen-india.in

