



# SW1470

## Control Software for NIM Power Supply Modules



## Features



- Java Program
- Supports Linux and Windows
- Up to 128 power supply channels easily managed
- USB or Ethernet (using the CAEN NIM8301 Crate) connection
- USB 2.0 Connection available on the N14xx front panel & Ethernet connection by the CAEN NIM8301 Crate Smart Fan Tray

## Description

**N14xx NIM Power Supply Modules** control can be performed remotely, via **USB** or **RS485**; the latter allows to build a N14XXs' daisy chain network, with the first module connected to the PC USB port and the subsequent ones daisy chained via RS485. It is therefore possible to build a daisy chain of up to **32 N14xx's (up to 128 channels)**.

The **CAEN NIM8301** Crate allows also to communicate with the first module in the chain also via Ethernet.

**SW1470 Control Software** is a Java Program supporting Linux and Windows that allows a complete and friendly management of the N14XX network.


### Drivers and Lib required

The USB connection (via Terminal emulator or via SW1470) requires driver of the N14xx USB chip controller (**FT232BM** - USB to serial UART interface) available at <http://www.ftdichip.com/FTDrivers.htm>. Since Java language does not provide a native support for the serial ports, a third-part library, the **RXTX library**, must be used for that purpose.

For more information about RXTX Library and USB driver click [here](#)

## Related Products

### N1470



4 Ch Reversible 8 kV/3 mA (8 W) NIM HV Power Supply Module (USB)

### N1471



4 Ch Reversible 5.5 kV/300  $\mu$ A NIM HV Power Supply Module (USB)

**C.A.E.N. SW1470 Control Software for NIM Power Supply Modules**

File Info About

C.A.E.N. N1471 4 Ch HV Power Supply V1.00 LocalBus Address: 0 Connected by:Ethernet IP: 10.0.7.111

Chan	VSet (V)	ISet (uA)	VMon (V)	IMon (uA)	Power	Status	Trip	Maxv (V)	RDWn (V/s)	RUp (V/s)	PWDn	Polarity	IMon Range
Chan 00	5000.0	310.000	0.0	0.000	Off	OFF	10.0	5600	500	500	KILL	-	HIGH
Chan 01	5500.0	310.000	0.0	0.000	Off	OFF	10.0	5600	500	500	KILL	-	HIGH
Chan 02	5500.0	310.000	0.0	0.000	Off	OFF	10.0	5600	500	500	KILL	-	HIGH
Chan 03	5500.0	310.000	0.0	0.000	Off	OFF	10.0	5600	500	500	KILL	-	HIGH

Internal Supply **OK** HV Clock [200KHz] Over Power **NO** Interlock Active **CLOSED** Local Bus Termination **OFF** Control Mode **REMOTE**

C.A.E.N. N1470A 2 Ch HV Power Supply V2.00 LocalBus Address: 1 Connected by:Ethernet IP: 10.0.7.111

Chan	VSet (V)	ISet (uA)	VMon (V)	IMon (uA)	Power	Status	Trip	Maxv (V)	RDWn (V/s)	RUp (V/s)	PWDn	Polarity	IMon Range
Chan 00	5000.0	3100.00	0.0	0.00	Off	KILL	10.0	8100	500	50	KILL	+	HIGH
Chan 01	5000.0	3100.00	0.0	0.00	Off	KILL	10.0	8100	500	50	KILL	+	HIGH

Internal Supply **OK** HV Clock [200KHz] Over Power **NO** Interlock Active **CLOSED** Local Bus Termination **OFF** Control Mode **REMOTE**

C.A.E.N. N1471 4 Ch HV Power Supply V1.00 LocalBus Address: 2 Connected by:Ethernet IP: 10.0.7.111

Chan	VSet (V)	ISet (uA)	VMon (V)	IMon (uA)	Power	Status	Trip	Maxv (V)	RDWn (V/s)	RUp (V/s)	PWDn	Polarity	IMon Range
Chan 00	5000.0	31.000											
Chan 01	5000.0	31.000											
Chan 02													
Chan 03													

Internal Supply **OK**

**SW1470 Info**  
[25/03/2010 10:50:58] Added new N1471 Board, LocalBus address 0, controlled by Ethernet IP: 10.0.7.111  
Communications OK:  
Added new N1470A Board, LocalBus address 1, controlled by Ethernet IP: 10.0.7.111  
Communications OK:  
Added new N1471 Board, LocalBus address 2, controlled by Ethernet IP: 10.0.7.111  
Communications OK:

**About SW1470**

**CAEN**  
Tools for Discovery

**SW1470** Version 1.0  
Written by: Alberto Lucchesi  
CAEN Computing Division  
A© C.A.E.N. SpA  
Info & support: support.computing@caen.it

Status: OK Communications

**This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit [www.caen.it](http://www.caen.it) © Caen S.p.A - 2024**

**CAEN S.p.A.**

Via Vetraia 11  
55049 - Viareggio  
Italy

**Phone +39.0584.388.398**

**Fax +39.0584.388.959**

**info@caen.it**

**www.caen.it**

**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt, Germany

**Phone +49 (0)212.254.4077**

**Mobile +49 (0)151.16.548.484**

**info@caen-de.com**

**www.caen-de.com**

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA

**Phone +1.718.981.0401**

**Fax +1.718.556.9185**

**info@caentechnologies.com**

**www.caentechnologies.com**

**CAENspa India Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri West  
Mumbai, Maharashtra 400053, India

**info@caen-india.in**

**www.caen-india.in**

