


A1422C/D

Low Noise Fast Rise Time Charge Sensitive Preamplifier



Features



- Fast, low noise inverting preamplifier
- Gain 5 mV/MeV(Si)
- Positive or negative input signals
- 8 input channels
- Up to 200 pF detector capacitance supported
- Up to 2000 V (positive or negative) detector bias voltage
- Low noise input stage composed JFETs diode protected

Description

The CAEN **Mod.A1422C/D** is a **5 mV/MeV** 8-channel preamplifier implemented into an alloy box and features **SHV** connectors for the DET/IN and HV/IN signals, **2 RJ45** for the differential outputs, **LEMO** connectors for the TEST IN and COM OUT and a cable with a **D-type 9** pin male connector for the power supply. The module accepts both positive and negative input charge pulses and provide a differential output of $\pm 5V$ max on 1 k Ω termination, and $\pm 2.5V$ on 100 Ω back termination. A test input accepts positive and negative signals for calibration purposes. Moreover the sum of the 8 energy outputs, with 1 μs decay time, is available on the COM OUT output.

A customization is available to improve amplifier protection from transients at detector input connection.

Technical Specifications

No. of Channels

A1422CDxxx: 8-channels

Packaging

Alloy Box; Dimensions WxHxD connector included: 100.0 x 50.0 x 162 mm³

Charge Sensitivity

5 mV/MeV (Si)

Recommended Range of Input Capacitance

F2: up to 200 pF

Output Linear Range

±5 V 1 kΩ termination ±2.5 V 100 Ω termination

COM Output Linear Range

±5 V 50 Ω termination

Open Loop Gain

5 mV/MeV (F2) $> 1 \cdot 10^5$

E²CRP Maximum energy-squared count-rate product

$E^2 \text{CRP}$: $1.57 \cdot 10^{10} \text{ MeV}^2/\text{s}$

Energy Sensitivity

5 mV/MeV (Si)

Decay Time

5 mV/MeV(Si) 50 μs

Integral non Linearity

$\leq \pm 0.05\%$ (0 ÷ ±10 V 1 kΩ termination)

Temperature Instability

$< \pm 100 \text{ ppm}/^\circ\text{C}$ (0 to 50°C)

Noise FWHM keV (Si)

Model Detector 0 pF Capacitance 200 pF F2 - 5 mV/MeV $< 4.7 < 7.6$

Rise Time

Model Detector 0pF Capacitance 200pF F2 - 5 mV/MeV $< 5 \text{ ns} < 15$

Detector Bias Voltage

± 2000 V max

Ordering Options

Code	Description
WA1422CD05F2	A1422CD05F2 - 8 Ch. Charge Preamplifier, 5mV/MeVgain, Cdet<200pF with differential output RoHS


Related Products

N5424



Quad NIM Power Distributor for A1422 Preamplicifiers

A1422H



Low-Noise Fast-Rise-Time Charge-Sensitive Preamplicifiers

V5425



Quad VME Power Distributor

DT5423



Quad Linear Desktop Power distributor for A1422 Preamplicifiers and DT57xx Digitizers

A1422



Low-Noise Fast-Rise-Time Charge-Sensitive Preamplicifiers (Boxed)

A1422E/F



Low Noise Fast Rise Time Charge Sensitive Preamplifiers

N1068




16 ch Programmable Spectroscopy Amplifier and 16ch CFD

N968



Spectroscopy Amplifier

Gallery



This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it © Caen S.p.A - 2024

CAEN S.p.A.

Via Vetraia 11
55049 - Viareggio
Italy

Phone +39.0584.388.398

Fax +39.0584.388.959

info@caen.it

www.caen.it

CAEN GmbH

Brunnenweg 9
64331 Weiterstadt, Germany

Phone +49 (0)212.254.4077

Mobile +49 (0)151.16.548.484

info@caen-de.com

www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA

Phone +1.718.981.0401

Fax +1.718.556.9185

info@caentechnologies.com

www.caentechnologies.com

CAENspa India Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri West
Mumbai, Maharashtra 400053, India

info@caen-india.in

www.caen-india.in

