



# A1480

## DC Power Input Equalizer for N14XX Family



## Features



- Optional DC Input Power Equalizer for N14XX Family

## Description

Optional DC Input Power Equalizer for **N14XX** Family: it allows to use also the NIM  $\pm 6V$  output to supply the units, increasing the total output power.

## Ordering Options

Code	Description	
WA1480XAAAAA	A1480 - DC Power Input Equalizer for N14XX Family	RoHS


## Related Products

### N1470ET



4 Ch Reversible 8 kV/3 mA (8 W) NIM HV Power Supply Module (USB/Ethernet/T.screen)

### N1470



4 Ch Reversible 8 kV/3 mA (8 W) NIM HV Power Supply Module (USB)

### N1471ET



4 Ch Reversible 5.5 kV/300  $\mu$ A NIM HV Power Supply Module (USB/Ethernet/T.screen)

### N1471H



4 Ch Reversible 5.5 kV/20  $\mu$ A NIM HV Power Supply High Accuracy Module (USB)

### N1471HET




4 Ch Reversible 5.5 kV/20  $\mu$ A NIM HV Power Supply High Accuracy Module (USB/Ethernet/T.screen)

## N1471



4 Ch Reversible 5.5 kV/300  $\mu$ A NIM HV Power Supply Module (USB)

---



**This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit [www.caen.it](http://www.caen.it) © Caen S.p.A - 2024**

**CAEN S.p.A.**

Via Vetraia 11  
55049 - Viareggio  
Italy

**Phone +39.0584.388.398**

**Fax +39.0584.388.959**

**info@caen.it**

**www.caen.it**

**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt, Germany

**Phone +49 (0)212.254.4077**

**Mobile +49 (0)151.16.548.484**

**info@caen-de.com**

**www.caen-de.com**

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA

**Phone +1.718.981.0401**

**Fax +1.718.556.9185**

**info@caentechnologies.com**

**www.caentechnologies.com**

**CAENspa India Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri West  
Mumbai, Maharashtra 400053, India

**info@caen-india.in**

**www.caen-india.in**

