


# A1481

## Kill Signal Adapter for N14XX series



## Description

The **A1481 Kill Signal Adapter** is a pluggable adapter that allows external DC level to enable the HV channel of the N14XX NIM HV Power Supply series.

The A1481 has on its rear side AMP 2 pin plugs to be inserted into the KILL inputs of the module's "Channel control section"; the front side has LEMO connectors.

This is necessary in order to use the modules with some detectors. Kill input is provided through LEMO connectors.

Two version are available, with a slightly different operation logic.

## Ordering Options

Code	Description	
WA1481XAAAAA	A1481 - N14xx Kill Signal Adapter	RoHS

## Related Products

### N1471



4 Ch Reversible 5.5 kV/300  $\mu$ A NIM HV Power Supply Module (USB)

### N1471H



4 Ch Reversible 5.5 kV/20  $\mu$ A NIM HV Power Supply High Accuracy Module (USB)

### N1470ET



4 Ch Reversible 8 kV/3 mA (8 W) NIM HV Power Supply Module (USB/Ethernet/T.screen)

### N1471ET



4 Ch Reversible 5.5 kV/300  $\mu$ A NIM HV Power Supply Module (USB/Ethernet/T.screen)

### N1470




4 Ch Reversible 8 kV/3 mA (8 W) NIM HV Power Supply Module (USB)

## N1471HET



4 Ch Reversible 5.5 kV/20  $\mu$ A NIM HV Power Supply High Accuracy Module  
(USB/Ethernet/T.screen)

---



**This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit [www.caen.it](http://www.caen.it) © Caen S.p.A - 2024**

**CAEN S.p.A.**

Via Vetraia 11  
55049 - Viareggio  
Italy

**Phone +39.0584.388.398**

**Fax +39.0584.388.959**

**info@caen.it**

**www.caen.it**

**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt, Germany

**Phone +49 (0)212.254.4077**

**Mobile +49 (0)151.16.548.484**

**info@caen-de.com**

**www.caen-de.com**

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA

**Phone +1.718.981.0401**

**Fax +1.718.556.9185**

**info@caentechnologies.com**

**www.caentechnologies.com**

**CAENspa India Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri West  
Mumbai, Maharashtra 400053, India

**info@caen-india.in**

**www.caen-india.in**

