


A315

Splitter



Description

The Mod. **A315** splits one input on two output signals. All the connectors are LEMO female type. The splitter is adapted for 50 Ω lines. The device is completely passive (no power supply is required); the amplitude on each output is one half of that on the input.

Ordering Options

Code	Description
WA315XAAAAAA	A315 - Splitter RoHS

Related Products

A313



20 dB Attenuator

DT5720A



A311



6 dB Attenuator

SP5607



Absorption Tool

A386



AC coupler BNC female to BCN male

A387



Input Filter cable BNC female to BNC male - 25 cm

A312




12 dB Attenuator

SP5600C



Educational Gamma Kit

Gallery



This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it © Caen S.p.A - 2024

CAEN S.p.A.

Via Vetraia 11
55049 - Viareggio
Italy

Phone +39.0584.388.398

Fax +39.0584.388.959

info@caen.it

www.caen.it

CAEN GmbH

Brunnenweg 9
64331 Weiterstadt, Germany

Phone +49 (0)212.254.4077

Mobile +49 (0)151.16.548.484

info@caen-de.com

www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA

Phone +1.718.981.0401

Fax +1.718.556.9185

info@caentechnologies.com

www.caentechnologies.com

CAENspa India Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri West
Mumbai, Maharashtra 400053, India

info@caen-india.in

www.caen-india.in

