

## **A318**

**Adapter for Clock  
signal FISCHER  
S101A004 male to  
3-pin AMPMODU IV  
female - 10 cm**



## Description

The adapter accepts NIM input and provides suitable output for the CLK-IN connector of the **CAEN Digitizer Series 1** and **741 Peak Sensing ADCs**.

## Ordering Options

Code	Description	
WA318XAAAAAA	A318 - Adapter for Clock signal FISCHER S101A004 male to 3-pin AMPMODU IV female - 10 cm	RoHS

## Related Products

### V1720



8 Input Channel 12bit 250 MS/s Digitizer

### DT5740D



32 Input Channel 12 bit 62.5MS/s Digitizer supporting DPP-QDC firmware

### V1740D



64 Input Channel 12 bit 62.5 MS/s Digitizer supporting DPP-QDC firmware

### V1724



8 Input Channel 14 bit 100 MS/s Digitizer

### VX1740



64 Input Channel 12bit 62.5 MS/s Digitizer

**N6724**



2/4 Channel 14 bit 100 MS/s Digitizer

**DT5740**



32 Input Channel 12 bit 62.5MS/s Digitizer

**DT5720**



4/2 Input Channel 12bit 250 MS/s Digitizer

**DT5724**



4/2 Input Channel 14 bit 100 MS/s Digitizer

**N6720**



4 Channel 12 bit 250 MS/s Digitizer

## N6740D



32 Channel 12bit 62.5 MS/s Digitizer

---

## N6740



32 Channel 12bit 62.5 MS/s Digitizer

---

## VX1740D



64 Input Channel 12bit 62.5 MS/s Digitizer supporting DPP-QDC firmware

---

## VX1720



8 Input Channel 12bit 250 MS/s Digitizer

---



**This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit [www.caen.it](http://www.caen.it) © Caen S.p.A - 2024**

**CAEN S.p.A.**

Via Vetraia 11  
55049 - Viareggio  
Italy

**Phone +39.0584.388.398**

**Fax +39.0584.388.959**

**info@caen.it**

**www.caen.it**

**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt, Germany

**Phone +49 (0)212.254.4077**

**Mobile +49 (0)151.16.548.484**

**info@caen-de.com**

**www.caen-de.com**

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA

**Phone +1.718.981.0401**

**Fax +1.718.556.9185**

**info@caentechnologies.com**

**www.caentechnologies.com**

**CAENspa India Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri West  
Mumbai, Maharashtra 400053, India

**info@caen-india.in**

**www.caen-india.in**

