


## A392

**Adapter 16x LEMO  
00 female (or MCX  
male) to 2.54mm  
34-pin female - 50  
cm / 1 m**



## Description

The **CAEN A392** is a passive adapter for analog signals from 16 LEMO connectors on coaxial cables, 50  $\Omega$  terminated (MCX version available) to header-type connector suitable for the inputs of a **V792**.

## Ordering Options

Code	Description	
7A392MCXAAAA	A392MCX - 16 Channel Cable Adapter (Flat to MCX) for V792, 50cm $\pm$ 10% cables	RoHS
7A392MCXMAAA	A392MCXM - 16 Channel Cable Adapter (Flat to MCX)for V792, 1m $\pm$ 10% cables	RoHS
WA392MCXAAAA	A392MCX - Adapter 16x MCX male to 2.54 mm 34-pin female - 50 cm	RoHS
WA392MCXMAAA	A392MCXM - Adapter 16x MCX male to 2.54 mm 34-pin female - 1 m	RoHS
WA392XAAAAAA	A392 - Adapter 16x LEMO 00 male to 2.54 mm 34-pin female - 50 cm	RoHS
WA392XMAAAAA	A392M - Adapter 16x LEMO 00 male to 2.54 mm 34-pin female - 1 m	RoHS

## Related Products


### V792



32 Channel Multievent QDC

---

### A371



Adapter 2.54mm 34x2 pin male to 1.27mm 68-pin ERNI SMC female

---



**This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit [www.caen.it](http://www.caen.it) © Caen S.p.A - 2024**

**CAEN S.p.A.**

Via Vetraia 11  
55049 - Viareggio  
Italy

**Phone +39.0584.388.398**

**Fax +39.0584.388.959**

**info@caen.it**

**www.caen.it**

**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt, Germany

**Phone +49 (0)212.254.4077**

**Mobile +49 (0)151.16.548.484**

**info@caen-de.com**

**www.caen-de.com**

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA

**Phone +1.718.981.0401**

**Fax +1.718.556.9185**

**info@caentechnologies.com**

**www.caentechnologies.com**

**CAENspa India Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri West  
Mumbai, Maharashtra 400053, India

**info@caen-india.in**

**www.caen-india.in**

