

A746E

**Patch panel 32x
LEMO 00 female to
1.27mm 68-pin ERNI
SMC female**



Description

The CAEN **Mod. A746E** is a passive desktop patch panel to adapt 32 LEMO connectors to the high-density ERNI connector matching the input of **DT5740, DT5740D**.

SMC 68 pin connector **LEMO connector**

Dual Row 68 pin connector for A746E LEMO type coaxial connector for A746E
Consult our **connectors reference page** for technical information.

Ordering Options

Code	Description	
WA746EXAAAAA	A746E - Patch panel 32x LEMO 00 female to 1.27 mm 68-pin ERNI SMC female	RoHS

Related Products

DT5740D



32 Input Channel 12 bit 62.5MS/s Digitizer supporting DPP-QDC firmware

DT5740



32 Input Channel 12 bit 62.5MS/s Digitizer

Gallery



This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it © Caen S.p.A - 2024

CAEN S.p.A.

Via Vetraia 11
55049 - Viareggio
Italy

Phone +39.0584.388.398

Fax +39.0584.388.959

info@caen.it

www.caen.it

CAEN GmbH

Brunnenweg 9
64331 Weiterstadt, Germany

Phone +49 (0)212.254.4077

Mobile +49 (0)151.16.548.484

info@caen-de.com

www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA

Phone +1.718.981.0401

Fax +1.718.556.9185

info@caentechnologies.com

www.caentechnologies.com

CAENspa India Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri West
Mumbai, Maharashtra 400053, India

info@caen-india.in

www.caen-india.in

