

## A7590D

**1 Ch. 80 V/50 mA  
Floating Digital  
Controlled SiPM  
Power Supply - PCB  
Mount**



## Features

- 2-75 V Output Voltage
- Programmable output ramp starting from 1.2V to output voltage
- 50 mA maximum output current
- 1 mV Output Step Voltage
- Low noise, 560  $\mu$ Vrms at 50mA
- Isolated output. It can be used to provide both positive and negative power supply
- I2C, UART, USB and Analog Voltage control
- HV SiPM temperature compensation with remote temperature sensing
- 100  $\mu$ V resolution voltage monitor
- 100 nA resolution current monitor
- 0.05°C resolution temperature monitor
- Programmable over-current protection
- Open source libraries in C/C++ and Python compatible with Arduino and Raspberry PI
- Working on Windows and Linux
- Zeus User interface to control multiple devices both via I2C, UART or USB
- THT and SMT compatible footprint

## Description

The **A7590D** is an isolated high voltage regulator specifically designed for SiPM bias. It has a built in temperature compensation controller with programmable coefficient.

The module can be digitally controlled by **UART** and **I2C**. Output voltage could be also controlled proportionally to an analog input, providing a voltage on a pin. A7590D can provide up to **50 mA** and the output voltage could be regulated between **20 V** and **80 V** with a resolution of **1 mV**. It is possible to set a programmable ramp starting from 1.2 V to the setpoint with configurable slope. The module is equipped with **24 bit ADC** to monitor output voltage and output current. A **PID** controller can be enabled to exploit the high resolution of the ADC to stabilize the output.

The module output is **electrically isolated** from the power supply and control interface. The module indeed can be used as floating power supply, or, grounding the negative or positive output, as a non-isolated positive or negative power supply.

## Technical Specifications

### Output Bias Voltage (Vset)

0 ÷ +2000 V DC

### Vset Resolution

1 V

### Vmon Resolution

1 V

### Ramp-Up/Ramp-Down

1 ÷ 500 V/s in steps of 1V

### Maximum Output Bias Current (Iset)

1 mA

## Ordering Options

Code	Description	
WA7590DUXAAA	A7590DU - USB Controlled Power Supply for SiPM 60V 70mA	RoHS
WA7590DXAAAA	A7590D - Digital Controlled Power Supply for SiPM 60V 70mA	RoHS

## Related Products

### DT5490



Desktop +85 V/10 mA Digital Controlled SiPM Power Supply with USB

---

**This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit [www.caen.it](http://www.caen.it) © Caen S.p.A - 2024**

**CAEN S.p.A.**

Via Vetraia 11  
55049 - Viareggio  
Italy

**Phone +39.0584.388.398**

**Fax +39.0584.388.959**

**info@caen.it**

**www.caen.it**

**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt, Germany

**Phone +49 (0)212.254.4077**

**Mobile +49 (0)151.16.548.484**

**info@caen-de.com**

**www.caen-de.com**

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA

**Phone +1.718.981.0401**

**Fax +1.718.556.9185**

**info@caentechnologies.com**

**www.caentechnologies.com**

**CAENspa India Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri West  
Mumbai, Maharashtra 400053, India

**info@caen-india.in**

**www.caen-india.in**

