

CAEN-CSA

Spectroscopy Amplifier Control Software



Features



- Innovative GUI
- Dashboard capability: allows to manage all the CAEN Spectroscopy Amplifiers
- Supports Linux and Windows
- Handles all the communication links provided by N568E and N1068 Spectroscopy Amplifiers:
 - Ethernet
 - USB

Description

Spectroscopy Amplifier Control Software is a software tool that allows an easy management of all the functional parameters of the CAEN Spectroscopy Amplifiers **N568E** and **N1068**. Module control can take place via either USB or TCP/IP and features an user friendly graphical user interface.

Supported operating systems are Windows and Linux.

Related Products

N568E






16 Channel Programmable Spectroscopy Amplifier (Low Noise)

N1068



16 ch Programmable Spectroscopy Amplifier and 16ch CFD



This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it © Caen S.p.A - 2024

CAEN S.p.A.

Via Vetraia 11
55049 - Viareggio
Italy

Phone +39.0584.388.398

Fax +39.0584.388.959

info@caen.it

www.caen.it

CAEN GmbH

Brunnenweg 9
64331 Weiterstadt, Germany

Phone +49 (0)212.254.4077

Mobile +49 (0)151.16.548.484

info@caen-de.com

www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA

Phone +1.718.981.0401

Fax +1.718.556.9185

info@caentechnologies.com

www.caentechnologies.com

CAENspa India Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri West
Mumbai, Maharashtra 400053, India

info@caen-india.in

www.caen-india.in

