



EASY_XML

EASY Rack Configuration Files



Features



- Board.xml configuration files for EASY Crate Configuration
- Contains operating parameters of All EASY boards
- CAEN EASY Rack Builder software required

Description

This page contains the .xml **Configuration Files** of **all EASY Boards**. These files contain the operating parameters of each EASY Board and are used to configure EASY Crates connected to the **A1676A** branch controller.

The A1676A branch controller is a single width board and will be housed in one SYx527 system, while the remote sections will operate in the "hostile area".

The User can build his crate configuration file using the **CAEN EASY Rack Builder**, a Java^(TM) application which allows EASY Users to create their customised EASY Crate configurations.

Related Products

A1676A



EASY 3000 Branch Controller

SY5527



Universal Multichannel Power Supply System / 19"wide, 4U-high (6 slot)

Gallery

The screenshot displays the 'EASY Rack Builder' software window. The title bar reads 'EASY Rack configuration About'. On the left, a list of 'Caen EASY Boards' includes models A3009, A3025, A3050, A3100, A3485, A3486, A3512N, A3512P, A3540, A3801, A3801A, A4601F, A4601H, and A4603. Below this list is a 'Selected Board' section with a vertical stack of buttons labeled A, 4, 6, 0, 3, where the '3' is highlighted in red. The main workspace shows six server racks arranged in a 2x3 grid, labeled 'EASY 3000 crate 01' through 'EASY 3000 crate 06'. Each rack has a 3D perspective view and a control panel with '48 VDC' and 'FAN CTRL' buttons. A blue arrow points from 'EASY4000' to 'EASY3000' between crates 01 and 02. At the bottom, a 'Branch Controller Name' field contains 'CAEN Branch Controller' with a note: '(require uppercase/lowercase letters and/or numbers.)'. A footer instruction reads 'Click on Board A4603 to drag it'.

This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it © Caen S.p.A - 2024

CAEN S.p.A.

Via Vetraia 11
55049 - Viareggio
Italy

Phone +39.0584.388.398

Fax +39.0584.388.959

info@caen.it

www.caen.it

CAEN GmbH

Brunnenweg 9
64331 Weiterstadt, Germany

Phone +49 (0)212.254.4077

Mobile +49 (0)151.16.548.484

info@caen-de.com

www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA

Phone +1.718.981.0401

Fax +1.718.556.9185

info@caentechnologies.com

www.caentechnologies.com

CAENspa India Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri West
Mumbai, Maharashtra 400053, India

info@caen-india.in

www.caen-india.in

