

NIM8306

2 Slot Switching 750 W Mini Crate



Features



- Cost effective, compact solution
- 12 A @ ± 6 V, 12A @ ± 12 V and 6 A @ ± 24 V output
- 2 air cooled NIM slots
- Ideal for small experiments, laboratories and universities
- Portable: 170 x 310 x 125 mm³ (WxHxD); 5 kg weight
- Equipped with long-life NIM connectors
- Wide AC input voltage range: 100 Vac to 240 Vac
- Active PFC
- Extra low noise switching power supply
- ON/OFF switch
- Power good and stand-by leds
- Short circuit protection

Description

The **CAEN Mod. NIM8306** is a 2 slot NIM crate with switching power supplies; the unit is provided with three ventilation fans.

The NIM8306 delivers 12A @ ± 6 V, 12A @ ± 12 V and 6A @ ± 24 V.

Safety features include: short circuit protection, Over / Under voltage protection, Over temperature protection. The Unit is powered by 100 \div 240 Vac, 50 \div 60 Hz, power factor 0.98% (230VAC).

Ordering Options

Code	Description
WNIM8306XAAA	NIM8306 - NIM crate, 2 slot, Switching PS ($\pm 6V$ 12A, $\pm 12V$ 12A, $\pm 24V$ 6A) RoHS

Related Products

NIM8305



2 Slot Switching 450 W Mini Crate



This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it © Caen S.p.A - 2024

CAEN S.p.A.

Via Vetraia 11
55049 - Viareggio
Italy

Phone +39.0584.388.398

Fax +39.0584.388.959

info@caen.it

www.caen.it

CAEN GmbH

Brunnenweg 9
64331 Weiterstadt, Germany

Phone +49 (0)212.254.4077

Mobile +49 (0)151.16.548.484

info@caen-de.com

www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA

Phone +1.718.981.0401

Fax +1.718.556.9185

info@caentechnologies.com

www.caentechnologies.com

CAENspa India Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri West
Mumbai, Maharashtra 400053, India

info@caen-india.in

www.caen-india.in

