

# SP5607

## Absorption Tool



## Features

- Spacers: one 4mm thick, five 10 mm thick
- Aluminium Absorbers: one 4mm thick, five 10 mm thick
- PMMA Absorbers: one 4 mm thick, five 10 mm thick

## Description

The Gamma Absorption Tool allows to perform gamma attenuation measurements. It is a modular tool and its design allows an easy connection to the **SP5606** bottom support. It is composed of:

- Spacers: one 4 mm thick, five 10 mm thick;
- Aluminium Absorbers: one 4 mm thick, five 10 mm thick;
- PMMA Absorbers: one 4 mm thick, five 10 mm thick.

## Ordering Options

Code	Description
WSP5607XAAA	SP5607 - Absorbition Tool <span data-bbox="1257 248 1337 300">RoHS</span>

## Related Products

### SP5609



Telescope Mechanics

---

### SP5600AN



Educational Kit - Premium Version

---

### SP5620CH



Cosmic Hunter

---

### SP5600D



Educational Beta Kit

---

**This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit [www.caen.it](http://www.caen.it) © Caen S.p.A - 2024**

**CAEN S.p.A.**

Via Vetraia 11  
55049 - Viareggio  
Italy

**Phone +39.0584.388.398**

**Fax +39.0584.388.959**

**info@caen.it**

**www.caen.it**

**CAEN GmbH**

Brunnenweg 9  
64331 Weiterstadt, Germany

**Phone +49 (0)212.254.4077**

**Mobile +49 (0)151.16.548.484**

**info@caen-de.com**

**www.caen-de.com**

**CAEN Technologies, Inc.**

1 Edgewater Street - Suite 101  
Staten Island, NY 10305  
USA

**Phone +1.718.981.0401**

**Fax +1.718.556.9185**

**info@caentechnologies.com**

**www.caentechnologies.com**

**CAENspa India Private Limited**

B205, BLDG42, B Wing,  
Azad Nagar Sangam CHS,  
Mhada Layout, Azad Nagar, Andheri West  
Mumbai, Maharashtra 400053, India

**info@caen-india.in**

**www.caen-india.in**

