


SP5609

**Telescope
Mechanics**



Features



Description

- Easy to mount

The **SP5609 Telescope Mechanics** is an additional tool for an easy setup of components that need to be on the same axis from different angles.

The SP5609 is composed by:

- n°1 rotary axis with desk support
- n°2 Clamp with screws
- n°2 angle bracket kit
- n°4 screws TPC 4X25 INOX CROSS DIN 7985
- n°4 zinc coated washers [4mm]
- n°4 zinc smooth washers
- n°4 inox nuts

Ordering Options

Code	Description
WSP5609XAAA	SP5609 - Telescope Mechanics RoHS

Related Products

SP5600AN



Educational Kit - Premium Version

SP5600D



Educational Beta Kit

SP5608




Scintillating Tile

SP5620CH



Cosmic Hunter

Gallery



This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it © Caen S.p.A - 2024

CAEN S.p.A.

Via Vetraia 11
55049 - Viareggio
Italy

Phone +39.0584.388.398

Fax +39.0584.388.959

info@caen.it

www.caen.it

CAEN GmbH

Brunnenweg 9
64331 Weiterstadt, Germany

Phone +49 (0)212.254.4077

Mobile +49 (0)151.16.548.484

info@caen-de.com

www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA

Phone +1.718.981.0401

Fax +1.718.556.9185

info@caentechnologies.com

www.caentechnologies.com

CAENspa India Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri West
Mumbai, Maharashtra 400053, India

info@caen-india.in

www.caen-india.in

