

SP5622

Detection system



Features



- Plastic scintillator (15 x 15 x 1cm²)
- Frontend electronic board
- AdvanSiD NUV-SiPM (4 x 4mm²) mounted in the tile corner at 45°

Description

The **SP5622** is the plastic scintillating tile for the Cosmic Hunter which support up to three tile for detection - coincidence of cosmic rays.

Ordering Options

Code	Description
WSP5622XAAA	SP5622 - Detection System RoHS


Related Products

SP5620CH



Cosmic Hunter

SP5622B



Portable scintillating tile for cosmic rays detection

Gallery



This document, or parts thereof, may not be reproduced in any form or by any means without written permission from Caen S.p.A. Although every effort has been made to ensure the accuracy of information presented in this catalog, Caen S.p.A reserves the right to modify its products specifications without giving any notice; for up to date information please visit www.caen.it © Caen S.p.A - 2024

CAEN S.p.A.

Via Vetraia 11
55049 - Viareggio
Italy

Phone +39.0584.388.398

Fax +39.0584.388.959

info@caen.it

www.caen.it

CAEN GmbH

Brunnenweg 9
64331 Weiterstadt, Germany

Phone +49 (0)212.254.4077

Mobile +49 (0)151.16.548.484

info@caen-de.com

www.caen-de.com

CAEN Technologies, Inc.

1 Edgewater Street - Suite 101
Staten Island, NY 10305
USA

Phone +1.718.981.0401

Fax +1.718.556.9185

info@caentechnologies.com

www.caentechnologies.com

CAENspa India Private Limited

B205, BLDG42, B Wing,
Azad Nagar Sangam CHS,
Mhada Layout, Azad Nagar, Andheri West
Mumbai, Maharashtra 400053, India

info@caen-india.in

www.caen-india.in

